

Stylos BA™


Bioactive Glass Matrix



Moldable.
Absorbant.
Irrigation Resistant.



Stylos BA™ | Bioactive Glass Matrix




Synergistic Effect

- An enhanced bone graft substitute composed of 45S5 bioactive glass, carbonate apatite anorganic bone mineral, and Type I collagen.
- A combination that drives a synergistic effect to induce osteoblasts proliferation, matrix maturation, and extracellular matrix mineralization.¹


Ideal Biologic Ensemble

- Uniform dispersion of bioactive glass and mineral particle throughout the collagen matrix mimics the characteristics of human bone, providing an ideal environment for bone regeneration.²



Bioactive Glass



- Cell Proliferation & Differentiation—45S5 Bioactive Glass stimulates growth and osteogenic differentiation of human primary osteoblasts.³
- 30% Bioactive Glass—formulated with the ideal amount for optimizing osteoblast growth.³
- 100-300µm Particle Size—provides a controlled rate of ion dissolution and surface reactivity for bone bonding and proliferation.^{4,5,6}




Similar Sized Macro & Micro Pores for Cell Migration

Carbonate Apatite

- Optimal Resorption & Remodeling—similar to human bone.^{7,8}
- Pores provide pathways for cell migration and attachment to lay down new bone.
- Higher osteoclastic and osteoblastic activity than BTCP and HA.⁹



Optimal Resorption Profile



Stylos BA™ | Bioactive Glass Matrix

Type I Collagen

- Highly absorbent, moldable, flexible, and resists migration upon irrigation.
- Binds proteins and cells and retains biological factors.¹⁰
- 100% resorbable through normal metabolic pathways.¹¹
- Intrinsic hemostatic properties control minor bleeding.^{11,12}



Stylos BA Strips

Part Number*	Part Description
STYS005	5cc, 3.2cm x 2cm x 0.8cm
STYS010	10cc, 6.25cm x 2cm x 0.8cm
STYS020	20cc, 12.5cm x 2cm x 0.8cm

*Additional Strip sizes available upon request



Stylos BA Putty

Part Number	Part Description
STYP002	2.5cc, 1.8cm x 1.8cm x 0.8cm
STYP005	5cc, 2.5cm x 2.5cm x 0.8cm
STYP010	10cc, 3.5cm x 3.5cm x 0.8cm

1. Ferreira, S.A., Young, G., Jones, J.R., Rankin, S. 2021. Bioglass/ Carbonate Apatite/ Collagen Composite Scaffold Dissolution Products Promote Human Osteoblast Differentiation. Materials Science & Engineering C, 118.
2. Data on file.
3. Xynos, I.D., Hukkanen, M.V., Batten, J.J., Buttery, L.D.K., Hench, L.L., Polak, J.M. (2000). Bioglass 45S5 stimulates osteoblast turnover and enhances bone formation In vitro: Implications and applications for bone tissue engineering. Calcif Tissue Int. 67(4), 321-9.
4. Oonishi, H., Kushitani, S., Yasukawa, E., Iwaki, H., Hench, L.L., Wilson, J., Tsuji, E., Sugihara, T. (1997). Particulate Bioglass Compared With Hydroxyapatite as a Bone Graft Substitute. Clinical Orthopaedics and Related Research, 334, 316-325, Lippincott-Raven Publishers, Philadelphia, PA.
5. Schepers, E.J.G., Ducheyne, P. (1997). Bioactive glass particles of narrow size range for the treatment of oral bone defects: a 1-24 month experiment.
6. Lindfors, N. C., Koski, I., Heikkilä, J. T., Mattila, K. and Aho, A. J. (2010), A prospective randomized 14-year follow-up study of bioactive glass and autogenous bone as bone graft substitutes in benign bone tumors. J. Biomed. Mater. Res., 94B, 157-164. doi:10.1002/jbm.b.31636
7. Matsuura, A., Kubo, T., Doi K., Hayashi, K., Morita, K., Yokota, R., Hayashi, H., Hirata, I., Okazaki, M., Akagawa, Y. (2009). Bone formation ability of carbonate.
8. Ellies, L.G., Carter, J.M., Natiella, J.R., Featherstone, J.D.B., Nelson, D.G.A. (1988). Quantitative analysis of early in vivo tissue response to synthetic apatite implants. J. of Biomed. Mater. Res., 22, 137-148.
9. Kanayama, K., Sriarj, W., Shimokawa, H., Ohya, K., Doi, Y., Shibutani, T. 2011. Osteoclast and Osteoblast Activities on Carbonate Apatite Plates in Cell Cultures. J. Biomaterials, 26, 435-436.
10. Geiger, M., Li, R.H., Friess, W. (2003). Collagen sponges for bone regeneration with rhBMP-2. Science Direct / Elsevier, 55, 1613-1629. <http://doi.org/10.1016/j.addr.2003.08.010>
11. Li, S.T. (2000). Biomedical Engineering Handbook, In JD Bronzino (Eds.), Biologic Biomaterials: Tissue Derived Biomaterials (Collagen) (1st ed.) 2, 42, 1-23, CRC Press, Boca Raton, FL.
12. Jaffe, R., Deykin, D. (1974). Evidence for a Structural Requirement for the Aggregation of Platelets by Collagen. The Journal of Clinical Investigation, 53, 875-883.

