

Orios[®]

Viable Bone Matrices



A Viable Autograft Alternative.



Orios® | Family of Viable Bone Matrices

The Orios family of viable bone matrices offers alternatives to autograft¹ that provide key elements ideal for bone formation:

- An osteoconductive three-dimensional scaffold with cortical and cancellous components.
- A demineralized bone scaffold with osteoinductive potential which provides exposure of signaling molecules and bone morphogenetic proteins.²
- Viable spine-derived cells for osteogenic supplementation.

Key Features & Benefits

- Allograft scaffold that provides a microenvironment for osteogenesis
- Viable cell population for osteogenic supplementation, with a minimum of 150,000 viable cells per cc of bone scaffold³
- Preparation with a proprietary DMSO-free cryoprotectant that allows for consistent delivery of viable allograft to the patient, with an average cell viability consistently exceeding 80% post-thaw³
- OR ease of use with no rinsing or decanting steps required, and 4-hour working window for implantation after cell thaw without loss of cell viability
- Availability in 3 unique scaffold blends to meet a wide variety of surgeon preferences and procedures: **Orios Moldable Plus**, **Orios Moldable**, **Orios Bone Matrix**.

Cell Preservation: A Different Approach

- Proper preservation of cellular allografts requires strict adherence to recovery and processing protocols.
- In each of the Orios viable bone matrices, viable spine-derived cells are collected from the vertebral body region of the donor and preserved with the use of a next-generation DMSO-free cryoprotectant.

Orios® | Family of Viable Bone Matrices



Figure 1:
DMSO-free cryoprotectant
coats the cells to prevent
crystalline damage




Figure 2:
Orios Representative
Sample





Figure 3:
2.5% DMSO Media
Sample

Evidence-based Application: Fusion

MIS-TLIF study demonstrated 96% fusion at 12 months.⁴



Operating Room Ease of Use

No rinsing or decanting
steps required

Average cell viability
consistently exceeds
80% post-thaw³

Minimum of
150,000 viable cells
per cc of allograft³

Four (4) hour working
window for implantation
after thaw without loss
of cell viability

Orios® | Family of Viable Bone Matrices



Orios Moldable Plus

Robust, fibrous moldability

Cortical shavings, crushed cancellous chips, and demineralized cortical bone microparticulate scaffold blend

Mixture creates a cohesive, fibrous consistency



Orios Moldable

For tight, defect packing

Cortical and cancellous bone microparticulate scaffold blend

Mixture creates a cohesive, wet sand consistency



Orios Bone Matrix

For moldable paste applications

Cortical and cancellous bone microparticulate scaffold blend with bone gel mixture

Mixture easily passes through a large or open bore syringe

Hydrophobic properties make it more resistant to lavage

| Part Number | Size |
|-------------|-------|
| BORMP025 | 2.5cc |
| BORMP005 | 5cc |
| BORMP010 | 10cc |
| BORMP015 | 15cc |

| Part Number | Size |
|-------------|-------|
| BORSM025 | 2.5cc |
| BORSM005 | 5cc |
| BORSM010 | 10cc |

| Part Number | Size |
|-------------|-------|
| BORS025 | 2.5cc |
| BORS005 | 5cc |
| BORS010 | 10cc |

1. MLR data on file at VIVEX Biologics, Inc.
2. Gruskin, E. et.al., Demineralized bone matrix in bone repair: history and use. Advanced Drug Delivery Reviews, 2012. 64:1063-1077.
3. Data on file at VIVEX Biologics, Inc.
4. Tally, William C, et al., Transforaminal Lumbar Interbody Fusion with Viable Allograft: 75 Consecutive Cases at 12-Month Follow-Up. International Journal of Spine Surgery, 2018. Vol. 12, No. 1 p. 76-84.

www.spinalelements.com

